

Parvalbumin alpha

Catalog No: #AP78910



Package Size: #AP78910-1 50ug #AP78910-2 100ug #AP78910-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Parvalbumin alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	2-110AA
Other Names	Pva
Accession No.	P02625Gene name:Pvalb
Uniprot	P02625
GeneID	25269;
Calculated MW	11.99
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

In muscle, parvalbumin is thought to be involved in relaxation after contraction. It binds two calcium ions.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only