

TYRO protein tyrosine kinase-binding protein

Catalog No: #AP78913



Package Size: #AP78913-1 50ug #AP78913-2 100ug #AP78913-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TYRO protein tyrosine kinase-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	22-113AA
Other Names	DAP12,KARAP,DNAX-activation protein 12,Killer-activating receptor-associated protein,KAR-associated protein
Accession No.	O43914Gene name:TYROBP
Uniprot	O43914
GenID	7305;
Calculated MW	10.12
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Non-covalently associates with activating receptors of the CD300 family. Cross-linking of CD300-TYROBP complexes results in cellular activation. Involved for instance in neutrophil activation mediated by integrin.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only