

# Pregnancy-associated glycoprotein 1

Catalog No: #AP78916



Package Size: #AP78916-1 50ug #AP78916-2 100ug #AP78916-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pregnancy-associated glycoprotein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Immunogen Description	198-297aa
Other Names	Pregnancy-specific protein B,PSP-B
Accession No.	Q29432Gene name:PAG1
Uniprot	Q29432
GeneID	281964;
Calculated MW	11
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Appears to be proteolytically inactive.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only