

## Transcobalamin-2

Catalog No: #AP78924



Package Size: #AP78924-1 50ug #AP78924-2 100ug #AP78924-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Transcobalamin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	Recombinant protein
Other Names	Transcobalamin II, TC II
Accession No.	Q9R0D6 Gene name: Tcn2
Uniprot	Q9R0D6
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Primary vitamin B12-binding and transport protein. Delivers cobalamin to cells.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only