

## GTPase HRas

Catalog No: #AP78931

Package Size: #AP78931-1 50ug #AP78931-2 100ug #AP78931-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GTPase HRas
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	Internal
Other Names	Hras1,H-Ras-1,Transforming protein p21,c-H-ras,p21ras
Accession No.	Q61411Gene name:Hras
Uniprot	Q61411
GeneID	15461;
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Ras proteins bind GDP/GTP and possess intrinsic GTPase activity.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only