

## Serine/threonine-protein kinase STK11

Catalog No: #AP78935



Package Size: #AP78935-1 50ug #AP78935-2 100ug #AP78935-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/threonine-protein kinase STK11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	Recombinant protein
Other Names	Lkb1,Liver kinase B1 homolog,LKB1
Accession No.	D4AE59Gene name:Stk11
Uniprot	D4AE59
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Isoform 2: Has a role in spermiogenesis.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only