

Thymic stromal lymphopoietin

Catalog No: #AP78940



Package Size: #AP78940-1 50ug #AP78940-2 100ug #AP78940-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thymic stromal lymphopoietin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	29-159AA
Other Names	null
Accession No.	Q969D9Gene name:TSLP
Uniprot	Q969D9
GeneID	85480;
Calculated MW	14.41
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Isoform 2: May act as an antimicrobial peptide in the oral cavity and on the skin.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only