

Prion-like protein doppel

Catalog No: #AP78941

Package Size: #AP78941-1 50ug #AP78941-2 100ug #AP78941-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prion-like protein doppel
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	27-152AA
Other Names	UNQ1830/PRO3443,DPL,PrPLP,Prion protein 2
Accession No.	Q9UKY0Gene name:PRND
Uniprot	Q9UKY0
GeneID	23627;
Calculated MW	13.86
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only