

## Actin, aortic smooth muscle

Catalog No: #AP78953

Package Size: #AP78953-1 50ug #AP78953-2 100ug #AP78953-3 1mg


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Actin, aortic smooth muscle
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	25-120AA
Other Names	GIG46,ACTSA,ACTVS,Alpha-actin-2,Cell growth-inhibiting gene 46 protein
Accession No.	P62736Gene name:ACTA2
Uniprot	P62736
GeneID	59;
Calculated MW	10.56
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Images



## Background

Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

---

Note: This product is for in vitro research use only