

Lipopolysaccharide-induced tumor necrosis factor-alpha factor homolog

Catalog No: #AP78961

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP78961-1 50ug #AP78961-2 100ug #AP78961-3 1mg

Description

Product Name	Lipopolysaccharide-induced tumor necrosis factor-alpha factor homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Immunogen Description	1-120AA
Other Names	SIMPLE, Small integral membrane protein of lysosome/late endosome
Accession No.	Q8QGW7 Gene name: LITAF
Uniprot	Q8QGW7
GeneID	374125;
Calculated MW	13.2
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Probable role in regulating transcription of specific genes.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only