

## Liver-expressed antimicrobial peptide 2

Catalog No: #AP78962



Package Size: #AP78962-1 50ug #AP78962-2 100ug #AP78962-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Liver-expressed antimicrobial peptide 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	37-76AA
Other Names	null
Accession No.	Q91V13Gene name:Leap2
Uniprot	Q91V13
GeneID	259301;
Calculated MW	4.4
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Has an antimicrobial activity.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only