

Interleukin-1 receptor-like 1

Catalog No: #AP78966

Package Size: #AP78966-1 50ug #AP78966-2 100ug #AP78966-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-1 receptor-like 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	19-328aa
Other Names	DER4,ST2,T1,Protein ST2
Accession No.	Q01638Gene name:IL1RL1
Uniprot	Q01638
GeneID	9173;
Calculated MW	34.1
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Isoform B: Inhibits IL-33 signaling.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only