

Cell surface A33 antigen

Catalog No: #AP78971



Package Size: #AP78971-1 50ug #AP78971-2 100ug #AP78971-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cell surface A33 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	22-235AA
Other Names	Glycoprotein A33
Accession No.	Q99795Gene name:GPA33
Uniprot	Q99795
GeneID	10223;
Calculated MW	23.54
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May play a role in cell-cell recognition and signaling.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only