

## Alpha-enolase

Catalog No: #AP79010

Package Size: #AP79010-1 50ug #AP79010-2 100ug #AP79010-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Alpha-enolase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	ENO1L1,MBPB1,MPB1,2-phospho-D-glycerate hydro-lyase,C-myc promoter-binding protein,Enolase 1,MBP-1,MPB-1,Non-neural enolase,Phosphopyruvate hydratase,Plasminogen-binding protein,NNE
Accession No.	P06733Gene name:ENO1
Uniprot	P06733
GeneID	2023;
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

MBP1 binds to the myc promoter and acts as a transcriptional repressor. May be a tumor suppressor.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only