

Visinin-like protein 1

Catalog No: #AP79011



Package Size: #AP79011-1 50ug #AP79011-2 100ug #AP79011-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Visinin-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	2-191AA
Other Names	Visl1,Neural visinin-like protein 1,NVL-1
Accession No.	P62761Gene name:Vsnl1
Uniprot	P62761
GeneID	26950;
Calculated MW	20.9
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Regulates (in vitro) the inhibition of rhodopsin phosphorylation in a calcium-dependent manner.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only