

ABI gene family member 3

Catalog No: #AP79013

Package Size: #AP79013-1 50ug #AP79013-2 100ug #AP79013-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABI gene family member 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-366aa
Other Names	NESH,New molecule including SH3,Nesh
Accession No.	Q9P2A4Gene name:ABI3
Uniprot	Q9P2A4
GeneID	51225;
Calculated MW	40.26
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May inhibit tumor metastasis (By similarity). In vitro, reduces cell motility.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only