

Matrix Gla protein

Catalog No: #AP79021



Package Size: #AP79021-1 50ug #AP79021-2 100ug #AP79021-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Matrix Gla protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	20-103AA
Other Names	null
Accession No.	P08494Gene name:Mgp
Uniprot	P08494
GeneID	25333;
Calculated MW	9.24
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Associates with the organic matrix of bone and cartilage. Thought to act as an inhibitor of bone formation.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only