

## Ropporin-1

Catalog No: #AP79030



Package Size: #AP79030-1 50ug #AP79030-2 100ug #AP79030-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ropporin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	C-terminal
Other Names	Rhophilin-associated 'oppo' protein,Rhophilin-associated protein 1
Accession No.	Q9ESG2Gene name:Ropn1
Uniprot	Q9ESG2
GeneID	76378;
Tag Info	His
Formulation	50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

null

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only