

CCDC102B Polyclonal Antibody

Catalog No: #40689

Package Size: #40689-1 50ul #40689-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CCDC102B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	CCDC102B Polyclonal Antibody detects endogenous levels of CCDC102B protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CCDC102B.
Target Name	CCDC102B
Other Names	CCDC102B; C18orf14; Coiled-coil domain-containing protein 102B
Accession No.	Swiss-Prot: Q68D86NCBI Gene ID: 79839
Uniprot	Q68D86
GeneID	79839;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using CCDC102B Polyclonal Antibody

Note: This product is for in vitro research use only