

Ornithine carbamoyltransferase, mitochondrial

Catalog No: #AP79033



Package Size: #AP79033-1 50ug #AP79033-2 100ug #AP79033-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ornithine carbamoyltransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	33-354AA
Other Names	Ornithine transcarbamylase,OTCase
Accession No.	P00480Gene name:OTC
Uniprot	P00480
GeneID	5009;
Calculated MW	35.42
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only