

Endogenous Bornavirus-like nucleoprotein 1

Catalog No: #AP79036

Package Size: #AP79036-1 50ug #AP79036-2 100ug #AP79036-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endogenous Bornavirus-like nucleoprotein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-360AA
Other Names	Endogenous Borna-like N element-1,EBLN-1
Accession No.	P0CF75Gene name:EBLN1
Uniprot	P0CF75
GeneID	340900;
Calculated MW	39.6
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May act as an RNA-binding protein. Highly homologous to the bornavirus nucleocapsid N protein that binds viral RNA and oligomerizes.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only