

N-acetylmuramoyl-L-alanine amidase

Catalog No: #AP79053



Package Size: #AP79053-1 50ug #AP79053-2 100ug #AP79053-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-acetylmuramoyl-L-alanine amidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant Protein
Other Names	UNQ3103/PRO10102,PGLYRPL,PGRPL,Peptidoglycan recognition protein 2,Peptidoglycan recognition protein long,PGRP-L
Accession No.	Q96PD5Gene name:PGLYRP2
Uniprot	Q96PD5
GeneID	114770;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May play a scavenger role by digesting biologically active peptidoglycan (PGN) into biologically inactive fragments. Has no direct bacteriolytic activity.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only