

Protein YIF1B

Catalog No: #AP79059

Package Size: #AP79059-1 50ug #AP79059-2 100ug #AP79059-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein YIF1B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-156AA
Other Names	PP4519, UNQ3073/PRO9905, YIP1-interacting factor homolog B
Accession No.	Q5BJH7 Gene name: YIF1B
Uniprot	Q5BJH7
GeneID	90522;
Calculated MW	17.16
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only