

# Neuronal-specific septin-3

Catalog No: #AP79069



Package Size: #AP79069-1 50ug #AP79069-2 100ug #AP79069-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Neuronal-specific septin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	58-331AA
Other Names	3-Sep
Accession No.	Q9UH03Gene name:42981
Uniprot	Q9UH03
GeneID	55964;
Calculated MW	30.14
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Filament-forming cytoskeletal GTPase (By similarity). May play a role in cytokinesis .

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only