

## H-2 class II histocompatibility antigen, A-B alpha chain

Catalog No: #AP79073



Package Size: #AP79073-1 50ug #AP79073-2 100ug #AP79073-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	H-2 class II histocompatibility antigen, A-B alpha chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	24-218AA
Other Names	null
Accession No.	P14434Gene name:H2-Aa
Uniprot	P14434
GeneID	14960;
Calculated MW	21.45
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

null

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only