

Multiple coagulation factor deficiency protein 2

Catalog No: #AP79075



Package Size: #AP79075-1 50ug #AP79075-2 100ug #AP79075-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Multiple coagulation factor deficiency protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	27-146AA
Other Names	SDNSF,Neural stem cell-derived neuronal survival protein
Accession No.	Q8NI22Gene name:MCFD2
Uniprot	Q8NI22
GeneID	90411;
Calculated MW	13.2
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

The MCFD2-LMAN1 complex forms a specific cargo receptor for the ER-to-Golgi transport of selected proteins. Plays a role in the secretion of coagulation factors.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only