

# Galectin-3-binding protein

Catalog No: #AP79079



Package Size: #AP79079-1 50ug #AP79079-2 100ug #AP79079-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Galectin-3-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	114-232aa
Other Names	M2BP,Basement membrane autoantigen p105,Lectin galactoside-binding soluble 3-binding protein,Mac-2-binding protein,Tumor-associated antigen 90K,MAC2BP
Accession No.	Q08380Gene name:LGALS3BP
Uniprot	Q08380
GeneID	3959;
Calculated MW	13.09
Tag Info	His
Formulation	50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Promotes intergrin-mediated cell adhesion. May stimulate host defense against viruses and tumor cells.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only