

Major facilitator superfamily domain-containing protein 2A

Catalog No: #AP79081

Package Size: #AP79081-1 50ug #AP79081-2 100ug #AP79081-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Major facilitator superfamily domain-containing protein 2A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	453aa-530aa
Other Names	HMFN0656,PP9177,UNQ300/PRO341,MFSD2
Accession No.	Q8NA29Gene name:MFSD2A
Uniprot	Q8NA29
GeneID	84879;
Calculated MW	8.58
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Play a role in thermogenesis via beta-adrenergic signaling pathway.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only