

Neuropeptide-like protein C4orf48 homolog

Catalog No: #AP79087



Package Size: #AP79087-1 50ug #AP79087-2 100ug #AP79087-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuropeptide-like protein C4orf48 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	29-90AA
Other Names	null
Accession No.	Q3UR78Gene name:Gm1673
Uniprot	Q3UR78
GeneID	381633;
Calculated MW	6.82
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only