

## Hepatoma-derived growth factor-related protein 3

Catalog No: #AP79102



Package Size: #AP79102-1 50ug #AP79102-2 100ug #AP79102-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hepatoma-derived growth factor-related protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	CGI-142,HDGF2,Hepatoma-derived growth factor 2,HDGF-2
Accession No.	Q9Y3E1Gene name:HDGFRP3
Uniprot	Q9Y3E1
GeneID	50810;
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Enhances DNA synthesis and may play a role in cell proliferation.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only