

Tsukushin

Catalog No: #AP79103



Package Size: #AP79103-1 50ug #AP79103-2 100ug #AP79103-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tsukushin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	UNQ850/PRO1788,E2IG4,LRRC54,TSK,E2-induced gene 4 protein,Leucine-rich repeat-containing protein 54
Accession No.	Q8WUA8Gene name:TSKU
Uniprot	Q8WUA8
GeneID	25987;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only