

Endogenous retrovirus group V member 1 Env polyprotein

Catalog No: #AP79106

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP79106-1 50ug #AP79106-2 100ug #AP79106-3 1mg

Description

Product Name	Endogenous retrovirus group V member 1 Env polyprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	ENVV1,HERV-V_19q13.41 provirus ancestral Env polyprotein 1
Accession No.	B6SEH8Gene name:ERVV-1
Uniprot	B6SEH8
GeneID	147664;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only