

Granulysin

Catalog No: #AP79110



Package Size: #AP79110-1 50ug #AP79110-2 100ug #AP79110-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Granulysin |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Human |
| Immunogen Description | Recombinant protein |
| Other Names | LAG2,NKG5,TLA519,Lymphokine LAG-2,Protein NKG5,T-cell activation protein 519 |
| Accession No. | P22749Gene name:GNLY |
| Uniprot | P22749 |
| GeneID | 10578; |
| Tag Info | His |
| Formulation | 50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0) |
| Storage | Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |

Background

Antimicrobial protein that kills intracellular pathogens. Active against a broad range of microbes, including Gram-positive and Gram-negative bacteria, fungi, and parasites. Kills Mycobacterium tuberculosis.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only