

Angiotensin-converting enzyme 2

Catalog No: #AP79111



Package Size: #AP79111-1 50ug #AP79111-2 100ug #AP79111-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	239-367AA
Other Names	UNQ5795/PRO19600,Angiotensin-converting enzyme 2
Accession No.	Q86XS5Gene name:ANGPTL5
Uniprot	Q86XS5
GeneID	253935;
Calculated MW	14.19
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only