

Asprosin

Catalog No: #AP79113



Package Size: #AP79113-1 50ug #AP79113-2 100ug #AP79113-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Asprosin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	2732-2871
Other Names	null
Accession No.	Gene name:ASPRO
Calculated MW	15.4
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Circulating asprosin, a protein hormone, responds to low dietary glucose by triggering the release of liver glucose stores, and the reduction of asprosin protects against the hyperinsulinism associated with metabolic syndrome.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only