

Geminin

Catalog No: #AP79118



Package Size: #AP79118-1 50ug #AP79118-2 100ug #AP79118-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Geminin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	CLTRNRLNRKRFRILG
Other Names	null
Accession No.	O75496Gene name:GMNN
Uniprot	O75496
GeneID	51053;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Inhibits the transcriptional activity of a subset of Hox proteins, enrolling them in cell proliferative control.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only