

Cytokine SCM-1 beta

Catalog No: #AP79120



Package Size: #AP79120-1 50ug #AP79120-2 100ug #AP79120-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Cytokine SCM-1 beta |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Human |
| Immunogen Description | Recombinant protein |
| Other Names | SCYC2,C motif chemokine 2,XC chemokine ligand 2 |
| Accession No. | Q9UBD3Gene name:XCL2 |
| Uniprot | Q9UBD3 |
| GeneID | 6846; |
| Tag Info | His |
| Formulation | 50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0) |
| Storage | Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |

Background

Chemotactic activity for lymphocytes but not for monocytes or neutrophils.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only