

Keratin, type II cytoskeletal 7

Catalog No: #AP79123



Package Size: #AP79123-1 50ug #AP79123-2 100ug #AP79123-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type II cytoskeletal 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	1-351aa
Other Names	Krt2-7,Cytokeratin-7,Keratin-7,Type-II keratin Kb7,CK-7,K7
Accession No.	Q9DCV7Gene name:Krt7
Uniprot	Q9DCV7
GeneID	110310;
Calculated MW	38.61
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Blocks interferon-dependent interphase and stimulates DNA synthesis in cells.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only