

High affinity immunoglobulin gamma Fc receptor I

Catalog No: #AP79139



Package Size: #AP79139-1 50ug #AP79139-2 100ug #AP79139-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	High affinity immunoglobulin gamma Fc receptor I
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	16-292AA
Other Names	FCG1,FCGR1,IGFR1,Fc-gamma RI,Fc-gamma RIA,FcRI,FcgammaRIa,CD64
Accession No.	P12314Gene name:FCGR1A
Uniprot	P12314
GeneID	2209;
Calculated MW	30.47
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

High affinity receptor for the Fc region of immunoglobulins gamma. Functions in both innate and adaptive immune responses.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only