

Fibronectin type III domain-containing protein 5

Catalog No: #AP79141


Package Size: #AP79141-1 50ug #AP79141-2 100ug #AP79141-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibronectin type III domain-containing protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-153aa
Other Names	FRCP2,Fibronectin type III repeat-containing protein 2
Accession No.	Q8NAU1Gene name:FNDC5
Uniprot	Q8NAU1
GeneID	252995;
Calculated MW	16.83
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

Irisin: Contrary to mouse, may not be involved in the beneficial effects of muscular exercise, nor in the induction of browning of human white adipose tissue.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only