

Cardiotrophin-like cytokine factor 1

Catalog No: #AP79151



Package Size: #AP79151-1 50ug #AP79151-2 100ug #AP79151-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cardiotrophin-like cytokine factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	BSF3,CLC,NNT1,B-cell-stimulating factor 3,Novel neurotrophin-1,BSF-3,NNT-1
Accession No.	Q9UBD9Gene name:CLCF1
Uniprot	Q9UBD9
GeneID	23529;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Cytokine with B-cell stimulating capability. Binds to and activates the ILST/gp130 receptor.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only