

Golgi membrane protein 1

Catalog No: #AP79157



Package Size: #AP79157-1 50ug #AP79157-2 100ug #AP79157-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Golgi membrane protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	40-205AA
Other Names	PSEC0242, UNQ686/PRO1326, C9orf155, GOLPH2, Golgi membrane protein GP73, Golgi phosphoprotein 2
Accession No.	Q8NBJ4 Gene name: GOLM1
Uniprot	Q8NBJ4
GeneID	51280;
Calculated MW	18.26
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Unknown. Cellular response protein to viral infection.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only