

Scavenger receptor cysteine-rich type 1 protein M130

Catalog No: #AP79170



Package Size: #AP79170-1 50ug #AP79170-2 100ug #AP79170-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Scavenger receptor cysteine-rich type 1 protein M130
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	46-401aa
Other Names	M130,Hemoglobin scavenger receptor,CD163
Accession No.	Q86VB7Gene name:CD163
Uniprot	Q86VB7
GeneID	9332;
Calculated MW	39.16
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

After shedding, the soluble form (sCD163) may play an anti-inflammatory role, and may be a valuable diagnostic parameter for monitoring macrophage activation in inflammatory conditions.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only