

T-cell surface glycoprotein CD5

Catalog No: #AP79192



Package Size: #AP79192-1 50ug #AP79192-2 100ug #AP79192-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T-cell surface glycoprotein CD5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	LEU1,Lymphocyte antigen T1/Leu-1,CD5
Accession No.	P06127Gene name:CD5
Uniprot	P06127
GeneID	921;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May act as a receptor in regulating T-cell proliferation.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only