

CD71 Polyclonal Antibody

Catalog No: #40705

Package Size: #40705-1 50ul #40705-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD71 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	CD71 Polyclonal Antibody detects endogenous levels of CD71 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human CD71 around the non-phosphorylation site of S24.
Target Name	CD71
Other Names	TFRC; Transferrin receptor protein 1; TR; TfR; TfR1; Trfr; T9; p90; CD antigen CD71
Accession No.	Swiss-Prot: P02786NCBI Gene ID: 7037
Uniprot	P02786
GeneID	7037;
SDS-PAGE MW	89kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using CD71 Polyclonal Antibody

Note: This product is for in vitro research use only