

Mucin-7

Catalog No: #AP79196



Package Size: #AP79196-1 50ug #AP79196-2 100ug #AP79196-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mucin-7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	CKLPPSPNNPPKFPNPHQPPKHPDK
Other Names	MG2,Apo-MG2,Salivary mucin-7
Accession No.	Q8TAX7Gene name:MUC7
Uniprot	Q8TAX7
GeneID	4589;
Tag Info	His
Formulation	50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May function in a protective capacity by promoting the clearance of bacteria in the oral cavity and aiding in mastication, speech, and swallowing. Binds *P.aeruginosa* pili.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only