

Versican core protein

Catalog No: #AP79201



Package Size: #AP79201-1 50ug #AP79201-2 100ug #AP79201-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Versican core protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	21-146AA
Other Names	CSPG2,Chondroitin sulfate proteoglycan core protein 2,Glial hyaluronate-binding protein,Large fibroblast proteoglycan,PG-M,Chondroitin sulfate proteoglycan 2,GHAP
Accession No.	P13611Gene name:VCAN
Uniprot	P13611
GeneID	1462;
Calculated MW	13.86
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May play a role in intercellular signaling and in connecting cells with the extracellular matrix. May take part in the regulation of cell motility, growth and differentiation. Binds hyaluronic acid.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only