

HLA class I histocompatibility antigen, B-27 alpha chain

Catalog No: #AP79212


Package Size: #AP79212-1 50ug #AP79212-2 100ug #AP79212-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HLA class I histocompatibility antigen, B-27 alpha chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	25-147AA
Other Names	HLAB,MHC class I antigen B*27
Accession No.	P03989Gene name:HLA-B
Uniprot	P03989
Calculated MW	13.53
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

Involved in the presentation of foreign antigens to the immune system.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only