

CD95 Polyclonal Antibody

Catalog No: #40709

Package Size: #40709-1 50ul #40709-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD95 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	CD95 Polyclonal Antibody detects endogenous levels of CD95 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human CD95.
Target Name	CD95
Other Names	FAS; APT1; FAS1; TNFRSF6; Tumor necrosis factor receptor superfamily member 6; Apo-1 antigen; Apoptosis-mediating surface antigen FAS; FASLG receptor; CD antigen CD95
Accession No.	Swiss-Prot: P25445NCBI Gene ID: 355
Uniprot	P25445
GeneID	355;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 823 cells using CD95 Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using CD95 Polyclonal Antibody

Note: This product is for in vitro research use only