

Angiotensin-related protein 2

Catalog No: #AP79234


Package Size: #AP79234-1 50ug #AP79234-2 100ug #AP79234-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin-related protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	180-287AA
Other Names	UNQ170/PRO196,ARP2,Angiotensin-like protein 2
Accession No.	Q9UKU9Gene name:ANGPTL2
Uniprot	Q9UKU9
GeneID	23452;
Calculated MW	11.88
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

Induces sprouting in endothelial cells through an autocrine and paracrine action.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only